


## Ursolic aldehyde

## Chemical Properties

CAS No. :	19132-81-1
Formula:	C <sub>30</sub> H <sub>48</sub> O <sub>2</sub>
Molecular Weight:	440.7
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Ursolic aldehyde can inhibit phosphatase of regenerating liver-3, with the IC(50) value of 50 +/- 0.3 mu M.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2691 mL	11.3456 mL	22.6912 mL
5 mM	0.4538 mL	2.2691 mL	4.5382 mL
10 mM	0.2269 mL	1.1346 mL	2.2691 mL
50 mM	0.0454 mL	0.2269 mL	0.4538 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481