


3-hydroxy-4-E-Hexenoic acid

Chemical Properties

CAS No. :	81357-32-6
Formula:	C ₆ H ₁₀ O ₃
Molecular Weight:	130.143
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-hydroxy-4-E-Hexenoic acid is a useful organic compound for research related to life sciences. The catalog number is T124181 and the CAS number is 81357-32-6.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.684 mL	38.4202 mL	76.8403 mL
5 mM	1.5368 mL	7.684 mL	15.3681 mL
10 mM	0.7684 mL	3.842 mL	7.684 mL
50 mM	0.1537 mL	0.7684 mL	1.5368 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481