

Polygalin J

Chemical Properties

CAS No. :

Formula: C41H46O21

Molecular Weight: 874.8

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Polygalin J can be used in research related to life sciences and its product number is TN6530.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1431 mL	5.7156 mL	11.4312 mL
5 mM	0.2286 mL	1.1431 mL	2.2862 mL
10 mM	0.1143 mL	0.5716 mL	1.1431 mL
50 mM	0.0229 mL	0.1143 mL	0.2286 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481