

## Carboxymethyl-PEG-Carboxymethyl (MW 5000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Carboxymethyl-PEG-Carboxymethyl (MW 5000) is a polyethylene glycol (PEG)-based linker designed specifically for the synthesis of PROTACs [1].
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands joined by a linker: one binds to an E3 ubiquitin ligase, and the other targets a specific protein. By leveraging the intracellular ubiquitin-proteasome system, PROTACs selectively degrade target proteins[1].

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481