


## Isopropylpiperazine

## Chemical Properties

CAS No. :	4318-42-7
Formula:	C7H16N2
Molecular Weight:	128.22
Storage:	Pure form: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Isopropylpiperazine, with CAS No. 4318-42-7, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Isopropylpiperazine provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.7991 mL	38.9955 mL	77.991 mL
5 mM	1.5598 mL	7.7991 mL	15.5982 mL
10 mM	0.7799 mL	3.8995 mL	7.7991 mL
50 mM	0.156 mL	0.7799 mL	1.5598 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481