

N,N'-DME-N,N'-Bis-PEG2-acid

Chemical Properties

CAS No. : 2062663-61-8


Formula: C18H36N2O8

Molecular Weight: 408.49

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	N,N'-DME-N,N'-Bis-PEG2-acid is a polyethylene glycol (PEG)-derived linker utilized in the synthesis of proteolysis targeting chimeras (PROTACs) [1].
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands connected by a linker: one ligand targets an E3 ubiquitin ligase, while the other targets the desired protein. They leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.448 mL	12.2402 mL	24.4804 mL
5 mM	0.4896 mL	2.448 mL	4.8961 mL
10 mM	0.2448 mL	1.224 mL	2.448 mL
50 mM	0.049 mL	0.2448 mL	0.4896 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481