


## Tetrahydro tanshinone I

## Chemical Properties

CAS No. :	126979-84-8
Formula:	C <sub>18</sub> H <sub>16</sub> O <sub>3</sub>
Molecular Weight:	280.32
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Tetrahydro tanshinone I is one compound of Fufang Danshen Prescription (FDP), liposome equilibrium dialysis system could simply and effectively estimate the absorption of multiple components in FDP in vivo, and predicts the potential bioactivity of these components, and thus contributes to the investigation of bioactive compounds with curative effect in FDP.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5674 mL	17.8368 mL	35.6735 mL
5 mM	0.7135 mL	3.5674 mL	7.1347 mL
10 mM	0.3567 mL	1.7837 mL	3.5674 mL
50 mM	0.0713 mL	0.3567 mL	0.7135 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481