

Bicyclo[1.1.1]pentane-1,3-dicarboxylic acid

Chemical Properties

CAS No. :	56842-95-6
Formula:	C7H8O4
Molecular Weight:	156.136
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Bicyclo[1.1.1]pentane-1,3-dicarboxylic acid ,with CAS No. 56842-95-6, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Bicyclo[1.1.1]pentane-1,3-dicarboxylic acid provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.4045 mL	32.0225 mL	64.0451 mL
5 mM	1.2809 mL	6.4045 mL	12.809 mL
10 mM	0.6405 mL	3.2023 mL	6.4045 mL
50 mM	0.1281 mL	0.6405 mL	1.2809 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481