


## 3-Methyl-5-phenylpyridazine

## Chemical Properties

CAS No. :	184021-10-1
Formula:	C <sub>11</sub> H <sub>10</sub> N <sub>2</sub>
Molecular Weight:	170.21
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-Methyl-5-phenylpyridazine, with CAS No. 184021-10-1, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 3-Methyl-5-phenylpyridazine provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.8751 mL	29.3755 mL	58.751 mL
5 mM	1.175 mL	5.8751 mL	11.7502 mL
10 mM	0.5875 mL	2.9375 mL	5.8751 mL
50 mM	0.1175 mL	0.5875 mL	1.175 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481