


Decanoic acid, barium salt

Chemical Properties

CAS No. :	13098-41-4
Formula:	C ₂₀ H ₃₈ BaO ₄
Molecular Weight:	479.84
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Decanoic acid, barium salt is a biochemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.084 mL	10.4201 mL	20.8403 mL
5 mM	0.4168 mL	2.084 mL	4.1681 mL
10 mM	0.2084 mL	1.042 mL	2.084 mL
50 mM	0.0417 mL	0.2084 mL	0.4168 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481