


Bay 65-1942 (R form)

Chemical Properties

CAS No. : 758683-21-5
 Formula: C₂₂H₂₅N₃O₄
 Molecular Weight: 395.45
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Bay 65-1942 (R form) is the less active enantiomer of Bay 65-1942, a selective and ATP-competitive IKK β inhibitor.
Targets(IC50)	Others,I κ B/IKK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5288 mL	12.6438 mL	25.2876 mL
5 mM	0.5058 mL	2.5288 mL	5.0575 mL
10 mM	0.2529 mL	1.2644 mL	2.5288 mL
50 mM	0.0506 mL	0.2529 mL	0.5058 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481