


N-Boc-diethanolamine

Chemical Properties

CAS No. :	103898-11-9
Formula:	C ₉ H ₁₉ NO ₄
Molecular Weight:	205.25
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	N-Boc-diethanolamine is a cleavable ADC linker used in the synthesis of antibody-drug conjugates (ADCs)[1].
Targets(IC50)	ADC Linker,PROTAC Linker
In vitro	ADCs are comprised of an antibody to which is attached an ADC cytotoxin through an ADC linker.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8721 mL	24.3605 mL	48.7211 mL
5 mM	0.9744 mL	4.8721 mL	9.7442 mL
10 mM	0.4872 mL	2.4361 mL	4.8721 mL
50 mM	0.0974 mL	0.4872 mL	0.9744 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Philip Stewart Low, et al. Small molecule ligand-targeted drug conjugates for anti-influenza chemotherapy and immunotherapy. W02020023323A1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481