

Trijugin C

Chemical Properties


CAS No. :

Formula: C₂₂H₃₄O₆

Molecular Weight: 394.508

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Trijugin C is a useful organic compound for research related to life sciences and the catalog number is T126176.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5348 mL	12.6739 mL	25.3479 mL
5 mM	0.507 mL	2.5348 mL	5.0696 mL
10 mM	0.2535 mL	1.2674 mL	2.5348 mL
50 mM	0.0507 mL	0.2535 mL	0.507 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481