

Compound T130777

Chemical Properties


CAS No. :

Formula: C₂₂H₃₄O₅

Molecular Weight: 378.51

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

1-[[5-(hydroxymethyl)-2,5,8a-trimethyl-1,4,4a,5,6,7,8,8a-octahydronaphthalen-1-yl]methyl]-4-methylidene-7-oxabicyclo[4.1.0]heptane-2,3,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6419 mL	13.2097 mL	26.4194 mL
5 mM	0.5284 mL	2.6419 mL	5.2839 mL
10 mM	0.2642 mL	1.321 mL	2.6419 mL
50 mM	0.0528 mL	0.2642 mL	0.5284 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481