

## Compound T129904

## Chemical Properties


CAS No. :

Formula: C48H72O14

Molecular Weight: 873.09

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(4R,4aR,5S,6aS,6bR,10S,12aR)-10-[(3,5-dihydroxy-4-[[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]oxan-2-yl)oxy]-5-hydroxy-4a-(hydroxymethyl)-2,2,6a,6b,9,9,12a-heptamethyl-1,2,3,4,4a,5,6,6a,6b,7,8,8a,9,10,11,12,12a,12b,13,14b-icosahydricen-4-yl benzoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1454 mL	5.7268 mL	11.4536 mL
5 mM	0.2291 mL	1.1454 mL	2.2907 mL
10 mM	0.1145 mL	0.5727 mL	1.1454 mL
50 mM	0.0229 mL	0.1145 mL	0.2291 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481