


Salverine

Chemical Properties

CAS No. : 6376-26-7
 Formula: C₁₉H₂₄N₂O₂
 Molecular Weight: 312.41
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Salverine is a bioactive chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2009 mL	16.0046 mL	32.0092 mL
5 mM	0.6402 mL	3.2009 mL	6.4018 mL
10 mM	0.3201 mL	1.6005 mL	3.2009 mL
50 mM	0.064 mL	0.3201 mL	0.6402 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481