

His Tag Immunoprecipitation Kit

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at 4°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description

TargetMol's His Tag Immunoprecipitation Beads specifically bind to His-tagged proteins and are suitable for immunoprecipitation (IP) of proteins, protein complexes, protein-nucleic acid complexes, and other antigens. The IP kit includes optimized, ready-to-use buffers that provide ideal reaction conditions for immunoprecipitation, enhancing the stability and consistency of the experiment. This product is suitable for antigens from cell lysates, cell culture supernatants, serum, ascitic fluid, and other biological samples.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481