

4-Hydrazinobenzenesulfonic acid-D4

Chemical Properties


CAS No. :

Formula: C₆D₄H₄N₂O₃S

Molecular Weight: 192.23

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	4-Hydrazinobenzenesulfonic acid-D4 is a deuterated compound of 4-Hydrazinobenzenesulfonic acid. 4-Hydrazinobenzenesulfonic acid has a CAS number of 98-71-5.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2021 mL	26.0105 mL	52.021 mL
5 mM	1.0404 mL	5.2021 mL	10.4042 mL
10 mM	0.5202 mL	2.6011 mL	5.2021 mL
50 mM	0.104 mL	0.5202 mL	1.0404 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481