


11R-CaN-AID

Chemical Properties

CAS No. :	650603-03-5
Formula:	C135H248N70O34S2
Molecular Weight:	3459.99
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	11R-CaN-AID is a potent, cell-permeant inhibitor of calcineurin [1] [2].
Targets(IC50)	Phosphatase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.289 mL	1.4451 mL	2.8902 mL
5 mM	0.0578 mL	0.289 mL	0.578 mL
10 mM	0.0289 mL	0.1445 mL	0.289 mL
50 mM	0.0058 mL	0.0289 mL	0.0578 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481