


2-(Pyren-1-ylaminocarbonyl)ethyl methanethiosulfonate

Chemical Properties

CAS No. :	384342-64-7
Formula:	C ₂₀ H ₁₇ NO ₃ S ₂
Molecular Weight:	383.48
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-(Pyren-1-ylaminocarbonyl)ethyl methanethiosulfonate serves as a fluorescent compound [1].
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6077 mL	13.0385 mL	26.077 mL
5 mM	0.5215 mL	2.6077 mL	5.2154 mL
10 mM	0.2608 mL	1.3038 mL	2.6077 mL
50 mM	0.0522 mL	0.2608 mL	0.5215 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481