


## 6-TET Dipivaloate

## Chemical Properties

CAS No. :	314734-87-7
Formula:	C <sub>31</sub> H <sub>24</sub> Cl <sub>4</sub> O <sub>9</sub>
Molecular Weight:	682.33
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	6-TET Dipivaloate is utilized as a nucleic acid sequencing probe.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4656 mL	7.3278 mL	14.6557 mL
5 mM	0.2931 mL	1.4656 mL	2.9311 mL
10 mM	0.1466 mL	0.7328 mL	1.4656 mL
50 mM	0.0293 mL	0.1466 mL	0.2931 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481