


## Acid yellow 199

## Chemical Properties

CAS No. :	70865-20-2
Formula:	C <sub>19</sub> H <sub>15</sub> N <sub>4</sub> NaO <sub>6</sub> S
Molecular Weight:	450.4
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Acid Yellow 199 is an azo dye that can be decolorized via a reduction reaction.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2202 mL	11.1012 mL	22.2025 mL
5 mM	0.444 mL	2.2202 mL	4.4405 mL
10 mM	0.222 mL	1.1101 mL	2.2202 mL
50 mM	0.0444 mL	0.222 mL	0.444 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481