


Ac-WEHD-pNA

Chemical Properties

CAS No. :	1135248-58-6
Formula:	C ₃₄ H ₃₇ N ₉ O ₁₁
Molecular Weight:	747.71
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ac-WEHD-pNA serves as a colorimetric substrate for caspases 1, 4, 5, and 14, widely used to assess caspase activity [1].
Targets(IC50)	Others,Caspase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3374 mL	6.6871 mL	13.3742 mL
5 mM	0.2675 mL	1.3374 mL	2.6748 mL
10 mM	0.1337 mL	0.6687 mL	1.3374 mL
50 mM	0.0267 mL	0.1337 mL	0.2675 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481