

ATTO 665

Chemical Properties

CAS No. : 1210944-09-4

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ATTO 665, a new fluorescent dye with high fluorescence yield, is suitable for most labeling applications with maximum excitation/emission wavelengths at 663/680 nm (in brackets).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481