

Avidin, chicken egg white

Chemical Properties

CAS No. : 1405-69-2

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Avidin, chicken egg white, is a glycoprotein obtained from egg protein with an exceptional affinity for biotin. This property makes it commonly used in combination with biotin for immunoassays to detect antigen locations in tissues [1].
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481