


Bamifylline

Chemical Properties

CAS No. :	2016-63-9
Formula:	C ₂₀ H ₂₇ N ₅ O ₃
Molecular Weight:	385.46
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Bamifylline, a drug in the xanthine chemical class, acts as a selective adenosine A1 receptor antagonist [1] [2].
Targets(IC50)	Adenosine Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5943 mL	12.9715 mL	25.943 mL
5 mM	0.5189 mL	2.5943 mL	5.1886 mL
10 mM	0.2594 mL	1.2972 mL	2.5943 mL
50 mM	0.0519 mL	0.2594 mL	0.5189 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481