


BCIP dipotassium

Chemical Properties

CAS No. :	390748-26-2
Formula:	C ₈ H ₄ BrClK ₂ NO ₄ P
Molecular Weight:	402.65
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	BCIP (dipotassium) is a versatile dye utilized in biological experiments to observe and analyze cell structures, track biomolecules, evaluate cell functions, distinguish cell types, detect biomolecules, study tissue pathology, and monitor microorganisms. Its applications extend from fundamental scientific research to a broad spectrum of clinical diagnostics, as well as traditional fields like textile dyeing and emerging areas such as functional textile processing, food pigments, and dye-sensitized solar cells.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4835 mL	12.4177 mL	24.8355 mL
5 mM	0.4967 mL	2.4835 mL	4.9671 mL
10 mM	0.2484 mL	1.2418 mL	2.4835 mL
50 mM	0.0497 mL	0.2484 mL	0.4967 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481