

Borocaptate sodium

Chemical Properties

CAS No. :	144885-51-8
Formula:	B12H12Na2S
Molecular Weight:	219.87
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Borocaptate sodium is an active molecule that can be used in life science related research. The CAS number of Borocaptate sodium is 144885-51-8.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5481 mL	22.7407 mL	45.4814 mL
5 mM	0.9096 mL	4.5481 mL	9.0963 mL
10 mM	0.4548 mL	2.2741 mL	4.5481 mL
50 mM	0.091 mL	0.4548 mL	0.9096 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481