


## BRD4 Inhibitor-30

## Chemical Properties

CAS No. :	1789731-20-9
Formula:	C <sub>28</sub> H <sub>38</sub> N <sub>6</sub> O <sub>4</sub>
Molecular Weight:	522.64
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	BRD4 Inhibitor-30 (Compound 1) acts as a BRD4 inhibitor, exhibiting an IC <sub>50</sub> value of 415 nM [1].
Targets(IC <sub>50</sub> )	Epigenetic Reader Domain

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9134 mL	9.5668 mL	19.1336 mL
5 mM	0.3827 mL	1.9134 mL	3.8267 mL
10 mM	0.1913 mL	0.9567 mL	1.9134 mL
50 mM	0.0383 mL	0.1913 mL	0.3827 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481