


C.I. Acid red 106

Chemical Properties

CAS No. :	6844-74-2
Formula:	C ₂₃ H ₁₇ N ₃ Na ₂ O ₉ S ₃
Molecular Weight:	621.57
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	C.I. Acid Red 106 is an acid red dye that is adsorbable and removable by bentone clay.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6088 mL	8.0441 mL	16.0883 mL
5 mM	0.3218 mL	1.6088 mL	3.2177 mL
10 mM	0.1609 mL	0.8044 mL	1.6088 mL
50 mM	0.0322 mL	0.1609 mL	0.3218 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481