


C.I. Acid yellow 49

Chemical Properties

CAS No. :	12239-15-5
Formula:	C ₁₆ H ₁₃ Cl ₂ N ₅ O ₃ S
Molecular Weight:	426.28
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	C.I. Acid Yellow 49 is an acidic yellow dye that can be adsorbed and removed by bentone [Bentone].
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3459 mL	11.7294 mL	23.4588 mL
5 mM	0.4692 mL	2.3459 mL	4.6918 mL
10 mM	0.2346 mL	1.1729 mL	2.3459 mL
50 mM	0.0469 mL	0.2346 mL	0.4692 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481