


Carbostyryl 165

Chemical Properties

CAS No. :	26078-23-9
Formula:	C ₁₂ H ₁₄ N ₂ O
Molecular Weight:	202.25
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Carbostyryl 165 is a neutral fluorescent dye (Ex: 351/364 nm, Em: >400 nm) [1].
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.9444 mL	24.7219 mL	49.4438 mL
5 mM	0.9889 mL	4.9444 mL	9.8888 mL
10 mM	0.4944 mL	2.4722 mL	4.9444 mL
50 mM	0.0989 mL	0.4944 mL	0.9889 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481