

CC-3240

## Chemical Properties

CAS No. : 2922675-91-8


Formula: C52H64N6O7S

Molecular Weight: 917.17

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	CC-3240 (compound 13), a molecular glue degrader of CaMKK2 based on CC-8977, has an IC50 of 9 nM [1].
Targets(IC50)	Molecular Glues,LRRK2

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0903 mL	5.4516 mL	10.9031 mL
5 mM	0.2181 mL	1.0903 mL	2.1806 mL
10 mM	0.109 mL	0.5452 mL	1.0903 mL
50 mM	0.0218 mL	0.109 mL	0.2181 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481