

COX-1/2-IN-6

Chemical Properties


CAS No. : 327024-84-0

Formula: C₂₂H₂₁NO₄

Molecular Weight: 363.41

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	COX-1/2-IN-6 (compound 4 h) acts as a dual inhibitor of COX-1 and COX-2, exhibiting potent inhibitory activity with IC ₅₀ values of 68 nM and 91 nM, respectively. It is utilized in the research of inflammatory diseases [1].
Targets(IC ₅₀)	COX

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7517 mL	13.7586 mL	27.5171 mL
5 mM	0.5503 mL	2.7517 mL	5.5034 mL
10 mM	0.2752 mL	1.3759 mL	2.7517 mL
50 mM	0.055 mL	0.2752 mL	0.5503 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481