


Cy3B

Chemical Properties

CAS No. :	228272-69-3
Formula:	C ₃₁ H ₃₂ N ₂ O ₆ S
Molecular Weight:	560.66
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Cy3B is an enhanced form of Cy3 (Cyanine3) dye, which is a fluorescent dye that emits light within the green to orange wavelength range.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7836 mL	8.9181 mL	17.8361 mL
5 mM	0.3567 mL	1.7836 mL	3.5672 mL
10 mM	0.1784 mL	0.8918 mL	1.7836 mL
50 mM	0.0357 mL	0.1784 mL	0.3567 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481