


Desacylсенegasaponin B

Chemical Properties

CAS No. : 163589-51-3
 Formula: C59H94O29
 Molecular Weight: 1267.36
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Desacylсенegasaponin B, an active compound extracted from the root of Polygala tenuifolia, exhibits potential anti-inflammatory activity [1].
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.789 mL	3.9452 mL	7.8904 mL
5 mM	0.1578 mL	0.789 mL	1.5781 mL
10 mM	0.0789 mL	0.3945 mL	0.789 mL
50 mM	0.0158 mL	0.0789 mL	0.1578 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481