


Disperse red 65

Chemical Properties

CAS No. :	16586-43-9
Formula:	C ₁₈ H ₁₈ ClN ₅ O ₂
Molecular Weight:	371.82
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Disperse Red 65, a red dye derived from anthraquinone, is utilized in plastics, textiles, and cosmetics.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6895 mL	13.4474 mL	26.8947 mL
5 mM	0.5379 mL	2.6895 mL	5.3789 mL
10 mM	0.2689 mL	1.3447 mL	2.6895 mL
50 mM	0.0538 mL	0.2689 mL	0.5379 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481