

Disperse yellow 54

Chemical Properties

CAS No. :	7576-65-0
Formula:	C ₁₈ H ₁₁ NO ₃
Molecular Weight:	289.28
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Disperse Yellow 54 is a multifunctional dye essential in biological experiments, aiding in the observation and analysis of cell structures, tracking of biomolecules, evaluation of cell functions, distinguishing cell types, detection of biomolecules, study of tissue pathology, and monitoring of microorganisms. Its applications span basic scientific research to clinical diagnostics, as well as traditional fields like textile dyeing and emerging areas such as functional textile processing, food pigments, and dye-sensitized solar cells.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4569 mL	17.2843 mL	34.5686 mL
5 mM	0.6914 mL	3.4569 mL	6.9137 mL
10 mM	0.3457 mL	1.7284 mL	3.4569 mL
50 mM	0.0691 mL	0.3457 mL	0.6914 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481