

E3 Ligase Ligand-linker Conjugate 28

Chemical Properties

CAS No. : 2229717-76-2


Formula: C₂₆H₃₃N₅O₆

Molecular Weight: 511.57

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	E3 Ligase Ligand-linker Conjugate 28, a combination of E3 ligase ligand and linker comprising Thalidomide and its corresponding linker, acts as a Cereblon ligand to recruit the CRBN protein and serves as a crucial intermediate in synthesizing complete PROTAC molecules.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9548 mL	9.7738 mL	19.5477 mL
5 mM	0.391 mL	1.9548 mL	3.9095 mL
10 mM	0.1955 mL	0.9774 mL	1.9548 mL
50 mM	0.0391 mL	0.1955 mL	0.391 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481