


GW759710A

Chemical Properties

CAS No. : 944795-23-7
 Formula: C₁₇H₁₄FN₅O
 Molecular Weight: 323.32
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	GW759710A is a FLuc inhibitor (IC ₅₀ : 0.73 μM) [1].
Targets(IC ₅₀)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0929 mL	15.4646 mL	30.9291 mL
5 mM	0.6186 mL	3.0929 mL	6.1858 mL
10 mM	0.3093 mL	1.5465 mL	3.0929 mL
50 mM	0.0619 mL	0.3093 mL	0.6186 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481