


Heteroclitin E

Chemical Properties

CAS No. : 140369-77-3
 Formula: C₂₇H₃₀O₉
 Molecular Weight: 498.52
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Heteroclitin E is a lignan derived from <i>Kadsura longipendunculata</i> [1].
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0059 mL	10.0297 mL	20.0594 mL
5 mM	0.4012 mL	2.0059 mL	4.0119 mL
10 mM	0.2006 mL	1.003 mL	2.0059 mL
50 mM	0.0401 mL	0.2006 mL	0.4012 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481