

## Hide powder Azure

## Chemical Properties

CAS No. : 37340-54-8

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Hide Powder Azure is a multifunctional dye. Used in biological experiments, it allows researchers to observe and analyze cell structures, track biomolecules, evaluate cell functions, distinguish cell types, detect biomolecules, study tissue pathology, and monitor microorganisms. Its applications span from basic scientific research to clinical diagnostics. Additionally, dyes play significant roles in traditional fields like textile dyeing and emerging areas such as functional textile processing, food pigments, and dye-sensitized solar cells.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481