


HPK1-IN-41

Chemical Properties

CAS No. : 2828455-23-6
 Formula: C₂₈H₃₃N₅O₂
 Molecular Weight: 471.59
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	HPK1-IN-41 (compound Z389) functions as an HPK1 inhibitor, exhibiting an IC ₅₀ value of 0.21 nM [1].
Targets(IC ₅₀)	MAPK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1205 mL	10.6024 mL	21.2049 mL
5 mM	0.4241 mL	2.1205 mL	4.241 mL
10 mM	0.212 mL	1.0602 mL	2.1205 mL
50 mM	0.0424 mL	0.212 mL	0.4241 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481