

HSD17B13-IN-34

Chemical Properties


CAS No. : 2770246-68-7

Formula: C₂₂H₁₃Cl₂F₃N₄O₃

Molecular Weight: 509.27

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	HSD17B13-IN-34 (75) functions as an inhibitor of 17β-Hydroxysteroid dehydrogenases (HSD17B13), exhibiting potency with an IC ₅₀ value of less than 0.1 μM against Estradiol. It is utilized primarily in the research of NAFLD (Nonalcoholic fatty liver diseases) [1].
Targets(IC ₅₀)	Dehydrogenase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9636 mL	9.818 mL	19.6359 mL
5 mM	0.3927 mL	1.9636 mL	3.9272 mL
10 mM	0.1964 mL	0.9818 mL	1.9636 mL
50 mM	0.0393 mL	0.1964 mL	0.3927 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481