


ICG Maleimide

Chemical Properties

CAS No. :	2143933-81-5
Formula:	C ₅₁ H ₅₆ N ₄ O ₆ S
Molecular Weight:	853.08
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	ICG Maleimide, a thiol-reactive near infrared (NIR) fluorescent dye, is used to generate a stable fluorescence signal in bioimaging.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1722 mL	5.8611 mL	11.7222 mL
5 mM	0.2344 mL	1.1722 mL	2.3444 mL
10 mM	0.1172 mL	0.5861 mL	1.1722 mL
50 mM	0.0234 mL	0.1172 mL	0.2344 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481