


Iopentol

Chemical Properties

CAS No. :	89797-00-2
Formula:	C ₂₀ H ₂₈ I ₃ N ₃ O ₉
Molecular Weight:	835.16
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Iopentol is a non-ionic contrast medium used for intravascular applications [1].
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1974 mL	5.9869 mL	11.9738 mL
5 mM	0.2395 mL	1.1974 mL	2.3948 mL
10 mM	0.1197 mL	0.5987 mL	1.1974 mL
50 mM	0.0239 mL	0.1197 mL	0.2395 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481