

Iso-isariin B

Chemical Properties


CAS No. : 2326436-93-3

Formula: C₃₀H₅₃N₅O₇

Molecular Weight: 595.77

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Iso-isariin B, an analogue of cyclocarboxyphenate, is utilized in antifungal studies [1].
Targets(IC50)	Antifungal

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6785 mL	8.3925 mL	16.785 mL
5 mM	0.3357 mL	1.6785 mL	3.357 mL
10 mM	0.1679 mL	0.8393 mL	1.6785 mL
50 mM	0.0336 mL	0.1679 mL	0.3357 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481