

KTt-45

Chemical Properties


CAS No. : 2244987-60-6

Formula: C₂₂H₂₄O₅

Molecular Weight: 368.42

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	KTt-45, a T-type calcium channel blocker, exhibits anticancer effects by inducing apoptosis in the HeLa cervical cancer cell line [1].
Targets(IC50)	Apoptosis, Calcium Channel, Caspase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7143 mL	13.5715 mL	27.1429 mL
5 mM	0.5429 mL	2.7143 mL	5.4286 mL
10 mM	0.2714 mL	1.3571 mL	2.7143 mL
50 mM	0.0543 mL	0.2714 mL	0.5429 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481