


Meteneprost potassium

Chemical Properties

CAS No. :	122576-55-0
Formula:	C ₂₃ H ₃₇ KO ₄
Molecular Weight:	416.64
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Meteneprost (U 46785) potassium acts as a cervical dilator by stimulating uterine contractions and widening the cervical canal.
Targets(IC50)	Others, Prostaglandin Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4002 mL	12.0008 mL	24.0015 mL
5 mM	0.480 mL	2.4002 mL	4.8003 mL
10 mM	0.240 mL	1.2001 mL	2.4002 mL
50 mM	0.048 mL	0.240 mL	0.480 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481