


p38- α MAPK-IN-6

Chemical Properties

CAS No. :	29368-40-9
Formula:	C ₁₅ H ₁₀ ClNO ₃
Molecular Weight:	287.70
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound 3a (p38- α MAPK-IN-6) is an inhibitor of the p38 α mitogen-activated protein kinase [1].
Targets(IC50)	p38 MAPK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4758 mL	17.3792 mL	34.7584 mL
5 mM	0.6952 mL	3.4758 mL	6.9517 mL
10 mM	0.3476 mL	1.7379 mL	3.4758 mL
50 mM	0.0695 mL	0.3476 mL	0.6952 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481